
Oxford[®] Comfort

Slings



Comfort

The Oxford® Comfort sling is an easy fit, contoured sling designed to accommodate the majority of clients. For use in conjunction with a 4-point positioning cradle, the patented Securi3 clip system ensures no inadvertent detachment of the sling.

- ✓ Available in polyester solid or net material options
- ✓ Optional padded leg sections for increased comfort
- ✓ Patented Securi3 clip attachment points
- ✓ Colour coded binding indicates sling size
- ✓ Positioning handles allow easy patient repositioning
- ✓ Modesty loop helps protect patient dignity during transfer
- ✓ Can be customised to fulfil patient specific needs

Size	Poly	Mesh	Padded	SWL
XS	LG6372	LG6382	LG6377	227 kg
S	LG6373	LG6383	LG6378	227 kg
M	LG6374	LG6384	LG6379	227 kg
L	LG6375	LG6385	LG6380	227 kg
XL	LG6376	LG6386	LG6381	227 kg



Chair to Chair



Bed to Chair/Chair to Bed



Floor to Bed/Chair



Chair to Commode/Toilet



Transfer to Bath



Rehab/Assisted Standing





At Enable we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



ENABLE LIFECARE PTY LTD

T 1300 370 370

enablelifecare.com.au