

Oxford[®] Deluxe Standing Slings




Deluxe Standing

The Oxford® Deluxe Standing sling is suitable for those patients who have a greater degree of weight bearing ability. It allows excellent access for toileting and is easy to fit, providing quick and effective transfers.


- ✓ Fully padded design offers excellent comfort
- ✓ Coloured attachment loops promote correct positioning
- ✓ Non-slip back section ensures the sling does not 'ride up'
- ✓ Adjustable waist strap for added support and security
- ✓ Can be customised to fulfil patient specific needs



Chair to Chair



Bed to Chair/Chair to Bed
(From/to seated position)



Chair to Commode/Toilet



Rehab/Assisted Standing



Size	Part No.	SWL
XS	SPECIAL	200 kg
S	SL13499N	200 kg
M	SL13500N	200 kg
L	SL13498N	200 kg
XL	SPECIAL	200 kg





At Enable we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



ENABLE LIFECARE PTY LTD

T 1300 370 370

enablelifecare.com.au