



## Empowering Australians through mobility and independence

ILS is a trusted Australian provider of assistive technology. Our highly trained and compassionate team are here to support you to make informed decisions.

We provide a wide range of equipment, with specialist teams who offer clinically led advice, ensuring equipment provides the best possible outcomes for clients.

## Quality Assurance

We are dedicated to providing the highest quality in assistive technology equipment, which is why all our Royale Medical products are registered with the Therapeutic Goods Administration (TGA) and certified to the Australian/New Zealand Standard (AS/NZS).

TGA Registration (ARTG: 291697)

AS/NZS 4688.2-2000 (R2016)

AS/NZS 4688.3-2000 (R2016)

AS/NZS 4688.4-2000 (R2016)

## Warranty

All Royale Medical Lift Chairs have a 2-year warranty on electrical components and 5-year structural warranty\*

\*Terms and conditions apply: 2-year electrical component warranty covers electrics including transformer, hand control, actuator, and all cables frame and steel lift/recline mechanism. 5-year structural warranty covers the timber frame.

Specifications	Mini	Standard
Overall Height:	102cm	110cm
Overall Width:	75cm	85cm
Overall Depth:	72cm	80cm
Armrest Height (from floor):	60.5cm	64cm
Seat Height (from floor):	44cm	50cm
Seat Width (measured on sitting area):	45cm	56cm
Seat Depth (measured on sitting area):	43cm	50cm
Safe Working Load (SWL):	150kg	150kg

# THE CHADWICK COLLECTION

Put yourself in complete control of your comfort with the best seat in the house.



PDT001-0422



# THE CHADWICK COLLECTION

## Innovation and Style

Exclusive to Independent Living Specialists, the Chadwick Collection features unique quad motor technology, with fully integrated, independent powered headrest, lumbar, backrest and footrest adjustments, putting you in complete control of your comfort and positioning needs.

Boasting a contemporary and elegant design, the Chadwick Collection suits all décor styles and was designed with comfort and ergonomics in mind.

## Clinical Benefits

With its powered headrest and lumbar support, the Chadwick Collection's postural adjustability can assist with accommodating an accentuated thoracic or lumbar kyphosis or lordosis.

Maintaining a normal curvature of the spine can help reduce the experience of back pain, as well as promote the ease of breathing and digestion.



## Key Features

- ✓ Soft vertical lift action to enable safe, easy movement and transfer out of the chair.
- ✓ Quad motor allowing for fully independent headrest, lumbar, backrest and footrest adjustments.
- ✓ 2-year warranty on all electrical components and 5-year structural warranty.\*
- ✓ Choice of 3 premium materials, each with colour options.
- ✓ Highest quality multi-layered foam for ultimate comfort and durability.
- ✓ Backup power supply for peace of mind.
- ✓ Hand control with unique home button bringing the chair back to its neutral position in one operation.
- ✓ USB port on the hand control to charge your electrical devices.

## Mini Range

Recommended for individuals 167cm/5'6" and below.

## Standard Range

Recommended for individuals 167cm/5'6" and above.

## Soft Touch Fabric

This scratch resistant material is breathable, warming and soft for comfort.



## Genuine Leather

This natural material has hypo-allergenic properties. It is durable, breathable, easy to maintain and attractive.



## Oxford Plush Leather

With all the benefits of Genuine Leather, this Italian leather is made with high quality hides producing a softer finish for increased comfort.

