

	The Select Range - In Studio Vinyl Available in Single Motor	The Signature Range - In Angora Fabric Available in Spacesaver Single Motor & Infinite Position Dual Motor		The Signature Range - In Angora Fabric Available in Infinite Position Dual Motor	The Luxury Range - Genuine Hide Leather Available in Spacesaver Single Motor & Infinite Position Dual Motor		The Luxury Range - Genuine Hide Leather Available in Infinite Position Dual Motor	The Chadwick Range - Soft Touch Fabric Available in Quad Motor	The Chadwick Range - Leather Available in Quad Motor	The Chadwick Range - Oxford Plush Leather Available in Quad Motor	The Chadwick Range - Soft Touch Fabric Available in Quad Motor	The Chadwick Range - Leather Available in Quad Motor	The Chadwick Range - Oxford Plush Leather Available in Quad Motor
ROYALE													
Motor	Select - Single Motor	Signature - Single Motor	Signature - Dual Motor	Signature - Dual Motor	Luxury - Single Motor	Luxury - Dual Motor	Luxury - Dual Motor	Infinite - Quad Motor	Infinite - Quad Motor	Infinite - Quad Motor	Infinite - Quad Motor	Infinite - Quad Motor	Infinite - Quad Motor
 Overall Height	105cm	105cm	105cm	98cm	105cm	105cm	98cm	110cm	110cm	110cm	102cm	102cm	102cm
 Overall Width	75cm	75cm	75cm	66cm	75cm	75cm	66cm	85cm	85cm	85cm	75cm	75cm	75cm
 Overall Depth	75cm	75cm	75cm	75cm	75cm	75cm	75cm	80cm	80cm	80cm	72cm	72cm	72cm
 Armrest Height (From floor)	62cm	62cm	62cm	56cm	62cm	62cm	56cm	64cm	64cm	64cm	60.5cm	60.5cm	60.5cm
 Seat Height (From floor)	48cm	48cm	48cm	45cm	48cm	48cm	45cm	50cm	50cm	50cm	44cm	44cm	44cm
 Seat Width (Measured on sitting area)	52cm	52cm	52cm	46cm	52cm	52cm	46cm	56cm	56cm	56cm	45cm	45cm	45cm
 Seat Depth (Measured on sitting area)	50cm	50cm	50cm	43cm	50cm	50cm	43cm	50cm	50cm	50cm	43cm	43cm	43cm
 Seat Max. Height (Lifted)	56cm	58cm	56cm	57cm	58cm	57cm	57cm	57cm	57cm	57cm	53.5cm	53.5cm	53.5cm
 Backrest Length (Fully reclined)	32cm	32cm	60cm	52cm	32cm	60cm	52cm	60cm	60cm	60cm	54cm	54cm	54cm
 Max. Legrest Height (Lifted)	46cm	46cm	48cm	47cm	46cm	48cm	47cm	52cm	52cm	52cm	47.5cm	47.5cm	47.5cm
 Safe Working Load (SWL)	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs	150kgs
Foam Density	28/70	28/70	28/70	28/70	28/70	28/70	28/70	35/70	35/70	35/70	35/70	35/70	35/70
Upholstery	Studio Vinyl	Premium Angora Fabric	Premium Angora Fabric	Premium Angora Fabric	Genuine Hide Leather Bovine Yearling	Genuine Hide Leather Bovine Yearling	Genuine Hide Leather Bovine Yearling	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)
Colour Options	 Charcoal black	 Symphony Blue	 Symphony Blue	 Symphony Blue	 Charcoal Black	 Charcoal Black	 Espresso Dark Brown	 Indigo Blue	 Taupe Beige	 Colonial Hemp Beige	 Indigo Blue	 Taupe Beige	 Jet Carbon Black
	 Vanilla Beige	 Dark Coffee Bean	 Dark Coffee Bean	 Dark Coffee Bean	 Espresso Brown	 Espresso Brown	 Taupe Beige	 Oyster Latte	 Espresso Dark Brown	 Jet Carbon Black	 Oyster Latte	 Espresso Dark Brown	 Marine Blue
	 Mushroom Brown	 Mushroom Grey	 Mushroom Grey		 Taupe Beige	 Taupe Beige		 Coffee Bean Brown		 Dark Cigar Brown	 Coffee Bean Brown		
	 Chocolate Brown									 Marine Blue			

**"BEST SEAT IN THE HOUSE"**