

Mayfair Collection

Chadwick Collection

Pressure Relieving Range



	Mayfair Collection							Chadwick Collection						Pressure Relieving Range		
Motor	Single	Single	Twin	Twin	Single	Twin	Twin	Quad	Quad	Quad	Quad	Quad	Quad	Twin	Twin	Twin
Material	Studio Vinyl	Premium Angora Fabric	Premium Angora Fabric	Premium Angora Fabric	Genuine Hide Leather Bovine Yearling	Genuine Hide Leather Bovine Yearling	Genuine Hide Leather Bovine Yearling	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)	Four Way Stretch Vapour Permeable PU Fabric	Four Way Stretch Vapour Permeable PU Fabric	Four Way Stretch Vapour Permeable PU Fabric
Size	Standard	Standard	Standard	Mini	Standard	Standard	Mini	Standard	Standard	Standard	Mini	Mini	Mini	Small	Medium	Large
Colours	<ul style="list-style-type: none"> ● Charcoal Black ● Vanilla Beige ● Mushroom Brown ● Chocolate Brown 	<ul style="list-style-type: none"> ● Symphony Blue ● Dark Coffee Bean ● Mushroom Grey 	<ul style="list-style-type: none"> ● Symphony Blue ● Dark Coffee Bean ● Mushroom Grey 	<ul style="list-style-type: none"> ● Symphony Blue ● Dark Coffee Bean 	<ul style="list-style-type: none"> ● Charcoal Black ● Espresso Brown ● Taupe Beige 	<ul style="list-style-type: none"> ● Charcoal Black ● Espresso Brown ● Taupe Beige 	<ul style="list-style-type: none"> ● Espresso Brown ● Taupe Beige 	<ul style="list-style-type: none"> ● Indigo Blue ● Oyster Latte ● Coffee Bean Brown 	<ul style="list-style-type: none"> ● Taupe Beige ● Espresso Brown 	<ul style="list-style-type: none"> ● Colonial Hemp Beige ● Jet Carbon Black ● Dark Cigar Brown ● Marine Blue 	<ul style="list-style-type: none"> ● Indigo Blue ● Coffee Bean Brown 	<ul style="list-style-type: none"> ● Taupe Beige ● Espresso Brown 	<ul style="list-style-type: none"> ● Jet Carbon Black ● Marine Blue 	<ul style="list-style-type: none"> ● Grey 	<ul style="list-style-type: none"> ● Grey 	<ul style="list-style-type: none"> ● Grey
Overall Height	1050mm	1050mm	1050mm	980mm	1050mm	1050mm	980mm	1100mm	1100mm	1100mm	1020mm	1020mm	1020mm	1060mm	1060mm	1060mm
Overall Width	750mm	750mm	750mm	660mm	750mm	750mm	660mm	850mm	850mm	850mm	750mm	750mm	750mm	750mm	750mm	820mm
Overall Depth	750mm	750mm	750mm	750mm	750mm	750mm	750mm	800mm	800mm	800mm	720mm	720mm	720mm	820mm	820mm	820mm
Armrest Height	620mm	620mm	620mm	560mm	620mm	620mm	560mm	640mm	640mm	640mm	605mm	605mm	605mm	620mm	640mm	640mm
Seat Height	480mm	480mm	480mm	450mm	480mm	480mm	450mm	500mm	500mm	500mm	440mm	440mm	440mm	480mm	520mm	520mm
Seat Width	520mm	520mm	520mm	460mm	520mm	520mm	460mm	560mm	560mm	560mm	450mm	450mm	450mm	500mm	500mm	595mm
Seat Depth	500mm	500mm	500mm	430mm	500mm	500mm	430mm	500mm	500mm	500mm	430mm	430mm	430mm	500mm	500mm	530mm
Backrest Length	320mm	320mm	600mm	520mm	320mm	600mm	520mm	600mm	600mm	600mm	540mm	540mm	540mm	450mm	450mm	450mm
Max. Legrest Height	460mm	460mm	480mm	470mm	460mm	480mm	470mm	520mm	520mm	520mm	475mm	475mm	475mm	500mm	500mm	500mm
Safe Working Load	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	160kg	160kg	160kg
Foam Density	28/70	28/70	28/70	28/70	28/70	28/70	28/70	35/70	35/70	35/70	35/70	35/70	35/70	35/70*	35/70*	35/70*

*Air bladder system to address pressure care needs.