PORTLAND

Dual Motor Riser Recliner Standard and Petite

> "Memory foam seat provides outstanding comfort"



The Portland Riser Recliner has been designed with comfort in mind. It features a memory-foam seat and the dual-motor mechanism enables the user to find their perfect comfort position.





PORTLAND


Dual Motor Riser Recliner Standard and Petite

> The Portland armchair features the luxurious Sofia fabric and is available in three colours: Russet, Mushroom and Graphite Grey.









- Memory foam seat upholstery provides added comfort and helps distribute pressure for long periods in chair.
- Three-tier waterfall back design provides extra comfort and support.
- Two motors enable the legrest to be elevated / lowered and the backrest to be reclined independently to give more flexible support
- Backrest can be positioned upright while the legrest is elevated; for example, to watch TV.
- It will rise and tilt users to their feet.
- Excellent snooze position for supreme comfort and relaxation.
- Quiet motors provide a smooth and discreet operation.
- Ergonomically contoured seat and back and padded upholstery for added comfort and support.
- Two convenient side pockets for easy storage.
- Safety key is built into hand controller to prevent accidental movement and entrapment.
- In-built fail safe battery back-up system which enables the chair to be raised and lowered in the event of a power failure*
- Comes with a simple to operate four button handset.
- Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility (e.g. those with MS, Arthritis, Rheumatism,
- Available in 3 colour options Russet, Mushroom and Graphite Grey.
- 2 year warranty**.
- Maximum weight capacity 150kg (23st/330lb).
 - * Batteries not included.

COLOUR OPTIONS

** Warranty conditions apply, for more information please see user manual.

Russet

PORTLAND

CLR22Y (Petite)

CODE: CLR12Y (Standard)

Max Weight Capacity: 150kg (23st/330lb)



Height: Std 110cm (44")



Width: Std 86cm (34") Pet 83cm (33")



Depth: Std 94cm (37") Pet 90cm (35")

Pet 101cm (40")



Seat Depth: Std 54cm (22") Pet 48cm (19")



Seat Height: Std 48cm (19") Pet 43cm (17")



Seat Width: Std 51cm (20") Pet 47cm (19")

www.drive-medical.com.au

Leaflet code: DDVH_010_Portland

Mushroom





Graphite Grey



