




# HERO+

## User Manual - S121

---

Adjustable  
Bed Cradle  
205kg



## Adjustable Bed Cradle Hero Medical 205kg



### Instructions for Assembly

1. Remove contents from packaging
2. Place receiver poles into base hole as shown in photos



## How to use the Adjustable Bed Cradle

This pole is designed to be used as an aid for bed transfer, please consult your healthcare professional for correct placement with different types of beds

## Adjusting the Height

Use pin button to raise or lower the bed cradle



**SWL : 205kg**

- This product is designed to be used daily in the home or community setting