

# HERO+


## User Manual - S12513

---

Dropside  
Bed Rail Tall  
160kg



Dropside Bed Rail Tall Hero Medical 160kg



## How to use the Dropside Bed Rail Tall

This product is designed to be used as a preventative aid for falling incurred from rolling off beds  
This product is designed for those with risk of falling from beds.

## Adjusting the Height

Adjust by using the black circular knob at the bottom of the rail, then adjust the height of the rail.



**SWL : 160kg**

- This product is designed to be used daily in the home or community setting