

Manual Wheelchairs

TiLite CR1




permobil




Key features

- Exceptionally lightweight frame
- ± 1 " Center of Gravity fine tuning adjustment
- Integrated, fold-down backrest for easy transport
- Individually made, just for you



Performance components

- Specially designed caster bearings for optimal maneuverability and weight savings
- PBO fiber and titanium reinforcements
- New carbon fiber caster wheels



Custom profiled sideguard

- Supports backrest without need for backrest hinges
- Individually designed for tailored fit




Dual-angled front end

- Designed to increase environmental access
- Optimized caster placement



[12.0096, 12.0161]





Component colors

Anodized options



COG bracket & trim rings



Specifications

| | |
|--------------------------|--------------|
| Transport weight | 9 lb |
| Seat width | 13" – 18" |
| Seat depth | 13" – 19" |
| Front seat height | 17.5" – 21" |
| Seat slope | 7°, 10°, 13° |
| Front frame angle | 85° |
| Seat back height | 10" – 17" |

| | |
|-----------------------------------|------------|
| Camber | 0°, 2°, 4° |
| Center of gravity | 2" – 5" |
| Backrest angle | 80° - 101° |
| Compatible with SmartDrive | Yes |
| Weight limit | 231 lb |
| Warranty | Limited |