Additional features found across the entire Zing product line such as leg abduction, swing-away tray, easy adjustments and comfortable transfer height further separates the Zing line from the competition, expands therapeutic opportunities and establishes a new standard of care for early intervention/pediatric standing.

PORTABLE, LIGHTWEIGHT AND FOLDABLE

- 3 foldable legs, a carrying handle and weight of under 20lbs, makes moving the stander from one room to another a breeze.
- With a small footprint when folded, traveling with your Zing stander is now possible and when a standing session is completed, the Zing Portable can be packed away without taking up much space.
- · Great for home, multi-user facility or school.

PLENTY OF GROWTH

- Fits children from infancy to 36" (91cm), and up to 36 lbs. (16kg).
- No minimum height or weight requirements which makes starting a standing program as young as possible a reality.
- Comfortable transfer height of 26" (66cm) with near flat loading makes transfers more manageable as children grow.

TRUE HIP PIVOT LEG ABDUCTION

- Pivot points are in line with the hips, promoting proper alignment in the hips, knees and feet.
- Feet maintain the angle with the rest of the legs no wedges needed and no additional adjustment to the knees or foot plates.
- Up to 15 degrees of abduction per leg accomplished through a single knob.
- · Promotes musculoskeletal development for children during critical growth years.
- Abduction can be adjusted and used throughout the entire positioning range.





EASY ADJUSTMENT

- Built in straps on the support pads.
- Tool-free adjustment of all components.
- Every adjustment can be made with a simple turn of a knob.

MULTI-USE TRAY-SUPINE

- The Supine features a multi-use tray that swings away either direction.
- Tray can remain attached during transfers.
- More activities are possible by swinging tray away during standing therapy.

NEW STANDARD OF CARE

Early intervention standing and the increasing prevalence of hip surveillance programs are leading to an evidence-based approach of standing children during the same window as typically developing peers. In other words, standing as young as possible for those children who are unable to do so unsupported is trending as a therapy intervention because of evidence being realized in many hip surveillance programs. Implementation of standing and standing in abduction is recommended in several hip surveillance programs as part of an intervention plan when hip integrity issues are identified. There are also a number of early intervention physical therapists advocating evidence-based approaches that include standing in abduction as a preventive approach to hip socket or acetabulum development prior to hip integrity diagnoses. Our bodies have evolved to stand early, stand often and children with developmental disabilities should not be exceptions in most cases.













The Zing Portable is lightweight (about 20 pounds) and features 3 fold-able legs that function like a traditional tripod. The ergonomic carry handle balances the weight of the unit and makes moving the stander from one place to another easy. Zing Portable was designed to fold up and pack away for either storage or travel with a minimal footprint.

zingstanders.com

Zing

Zing specifications

Portable	
Height Range	Up to 36" (91cm)
Weight Limit	Up to 36 lbs. (16.3kg)
Weight of basic unit (without tray or mast top)	14 lbs. (6.4kg)
Weight of Supine Mast Top w/tray and headrest	5.5 lbs. (2.5kg)
Weight of Prone Mast Top w/tray	2 lbs. (1 kg)
Tripod Frame Footprint, approximately	39"W x 32"L (91x81cm)
Transfer Height (Horizontal at hip pivot)	26" (66cm)
Leg Abduction	0 to 15°
Foot Plate Size	6"x2.75" (15x6cm)
Foot Height Range (Center of Pelvic Support to Footplate)	8.5"-17.5" (22-44cm)
Kneepad Inside Height, Width is variable (form-to-fit)	4" (10cm)
Kneepads Lateral Adjustment, per side	1.25" (3cm)
Kneepads Height Range, (Center of pelvic pad to center of knee)	
	4.5"-13.5" (11-34cm)
Foot straps- length over arch w/min 1" Velcro overlap	8.5" (21cm)
Tray Size Molded	7.5"x 12.5" w/3"Elbow Extensions (19 x 32 cm w/8cm Elbow Extensions)
Pelvic, Upper Body, and Head Support Pad Sizes	3.5"x8" (9x20cm)



AltimateMedical 262 W. 1st Street Morton, MN 56270



800.342.8968 zingstanders.com info@zingstanders.com

Zing Portable

zingstanders.com



Standing made possible anyplace, anytime.

Building upon the industry-leading innovation and features of the Zing multi-position stander, the new Zing Portable stander is a leap forward in product design. This new early intervention stander is the first "right-sized" product for standing children as young as possible (fits children from infancy to 36" and up to 36 pounds). Typically developing children will usually begin the processes of learning to stand between 9-13 months. The Zing Portable is designed to fit this age group better than any other product available, but that is just one of a long list of features that creates a new standard of care.

PORTABLE

LIGHTWEIGHT

FOLDABLE