

# Koala

Ready or not...



permobil



Dynamite performance in a small package means you can safely go wherever the adventure takes you.



### Miniflex Power Seating

Unique power seat elevation and tilt designed to provide valuable positioning and exceptional low seat to floor height.



### Parental Approved

Optional features such as front bumper and remote stop give any parent that extra piece of mind for their little explorer.

### Colour Options



Volcano Red



Popstar Pink



Extreme Yellow



Green Apple



Cobalt Blue



Purple Rain



Platinum Silver



Onyx Black

### Additional Hub Colour Options



Chili Pepper Red



Titan Blue



Phantom Gray



Vintage Wine



Crystal Blue



Sublime Green

### Specifications

Maximum User Weight	35 kg
Maximum Speed	6.5 km/h
Range	24 km*
Base Width	23"
Base Length	32"
Minimum Turning Radius	21"
Ground Clearance	2.5"
Recommended Battery Type	M40
Weight Including Batteries	102.5 kg
Drive Electronics	R-Net
Maximum Slope	10° up, 6° down†
Seat-to-Floor Height	15"

Optional Seat Elevator	8"
Power Tilt	0° – 25°
Manual Recline	85° – 120°
Manual Legrest Elevation	90° – 150°
Backrest Height	14", 16", 20"
Seat Depth	10", 12", 14"
Seat Width	10", 12", 14"
Armrest Pad Lengths	8", 10"

\*Actual driving range will vary based on driving and battery conditions.  
†Refer to Owner Manual for safe driving practices.