

Mayfair Collection

Chadwick Collection

Pressure Relieving Range



Size	Standard	Standard	Mini	Standard	Mini	Standard	Standard	Standard	Mini	Mini	Mini	Small	Medium	Large
Material	Studio Vinyl	Premium Angora Fabric	Premium Angora Fabric	Genuine Hide Leather Bovine Yearling	Genuine Hide Leather Bovine Yearling	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)	Soft Touch Fabric	Genuine Leather	Oxford Plush Leather (Italian)	Four Way Stretch Vapour Permeable PU Fabric	Four Way Stretch Vapour Permeable PU Fabric	Four Way Stretch Vapour Permeable PU Fabric
Model	Mayfair Select	Mayfair Signature	Mayfair Signature	Mayfair Luxury	Mayfair Luxury	Chadwick	Chadwick	Chadwick	Chadwick	Chadwick	Chadwick	Air Pressure Relieving	Air Pressure Relieving	Air Pressure Relieving
Motor	Single	Twin	Twin	Twin	Twin	Quad	Quad	Quad	Quad	Quad	Quad	Twin	Twin	Twin
Colours	<ul style="list-style-type: none"> Charcoal Black Vanilla Beige Mushroom Brown Chocolate Brown 	<ul style="list-style-type: none"> Symphony Blue Dark Coffee Bean Mushroom Grey 	<ul style="list-style-type: none"> Symphony Blue Dark Coffee Bean Taupe Beige 	<ul style="list-style-type: none"> Charcoal Black Espresso Brown Taupe Beige 	<ul style="list-style-type: none"> Espresso Brown Taupe Beige 	<ul style="list-style-type: none"> Indigo Blue Oyster Latte Coffee Bean Brown 	<ul style="list-style-type: none"> Taupe Beige Espresso Brown 	<ul style="list-style-type: none"> Colonial Hemp Beige Jet Carbon Black Dark Cigar Brown Marine Blue 	<ul style="list-style-type: none"> Indigo Blue Coffee Bean Brown 	<ul style="list-style-type: none"> Taupe Beige Espresso Brown 	<ul style="list-style-type: none"> Jet Carbon Black Marine Blue 	<ul style="list-style-type: none"> Grey 	<ul style="list-style-type: none"> Grey 	<ul style="list-style-type: none"> Grey
Overall Height	105cm	105cm	98cm	105cm	98cm	110cm	110cm	110cm	102cm	102cm	102cm	106cm	106cm	106cm
Overall Width	85cm	85cm	77cm	85cm	77cm	85cm	85cm	85cm	75cm	75cm	75cm	75cm	75cm	82cm
Overall Depth	85cm	85cm	80cm	85cm	80cm	80cm	80cm	80cm	72cm	72cm	72cm	82cm	82cm	82cm
Armrest Height	62cm	62cm	56cm	62cm	56cm	64cm	64cm	64cm	60.5cm	60.5cm	60.5cm	62cm	64cm	64cm
Seat Height	48cm	48cm	45cm	48cm	45cm	50cm	50cm	50cm	44cm	44cm	44cm	48cm	52cm	52cm
Seat Width	52cm	52cm	47cm	52cm	47cm	56cm	56cm	56cm	45cm	45cm	45cm	50cm	50cm	59.5cm
Seat Depth	46cm	50cm	43cm	50cm	43cm	50cm	50cm	50cm	43cm	43cm	43cm	50cm	50cm	53cm
Max. Legrest Height	46cm	48cm	47cm	48cm	47cm	52cm	52cm	52cm	47.5cm	47.5cm	47.5cm	50cm	50cm	50cm
Safe Working Load	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	150kg	160kg	160kg	160kg
Foam Density	28/70	28/70	28/70	28/70	28/70	35/70	35/70	35/70	35/70	35/70	35/70	35/70*	35/70*	35/70*

*Air bladder system to address pressure care needs.



Lift Chair Comparison Chart

