

# ROYALE

## Quick Shift Swivel

**Sit-to-Stand Transfer Aid**

Efficient support for safe and comfortable patient transfers.



A sit-to-stand transfer aid for users who can bear weight and push through their legs. Enables safe, comfortable transfers while reducing caregiver strain. Perfect for both home and professional care settings.

# Quick Shift Swivel

## Why Choose the Royale Quick Shift Swivel?

- **Supportive and Stable:** Sturdy steel frame provides essential stability during transfers, minimising risk and strain for users and caregivers – ensuring smooth, supported transfers.
- **Ergonomic and Adjustable:** The upper handle and lower leg supports are fully adjustable, promoting correct posture and customisable comfort.
- **Easy Foot Positioning:** The base plate includes clear left/right markers for precise foot placement, with a counter-balance support for caregivers.
- **Compact and Portable:** Durable steel frame with rear castors for easy transport and compact storage.
- **Versatile in Any Environment:** Suitable for homes or facilities, with an optional standing belt for added security.
- **Hygienic and Easy to Clean:** Wipe-down surfaces make it practical for daily use in both care facilities and home settings.



Padded Lower Leg Supports



Low Profile Base Plate  
Counterbalance Position for Attende

## Specifications

	S332021
Overall Width	41cm
Overall Depth	45cm
Overall Height	124cm
Height of Base Plate	2cm
Base Plate Width	41cm
Base Plate Depth	43cm
Padded Leg Support (W x H)	16cm x 19cm
Padded Leg Support Height Adjustment	34cm - 55cm
Handle Height (from base plate - bottom/top)	109cm - 125cm
Unit Weight	21kg
Castor Diameter	10cm
Safe Working Load (SWL)	150kg