

ROYALE

Quick Shift Deluxe

Sit-to-Stand Transfer Aid

Premium transfer aid with customisable features for comfort and control.



A flexible sit-to-stand transfer aid that supports active standing and short-distance transfers between seats. Enables caregivers to safely guide users between beds, chairs, and other seated positions.

Quick Shift Deluxe

Why Choose the Royale Quick Shift Deluxe?

- **Advanced Support System:** Multiple ergonomic grip points and a versatile upper handle bar provide enhanced stability for both user and caregiver.
- **Fully Customisable Comfort:** Tool-less adjustments for leg support height, width, and angle, with soft padding for user comfort.
- **Intelligent Base Design:** Low-profile base with recessed front fits between wheelchair wheels while maintaining a spacious foot platform.
- **Enhanced Safety Features:** Central locking system with colour-coded brake pedals for secure positioning.
- **Professional-Grade Durability:** Robust steel frame designed for daily use in clinical and home environments.
- **Complete Transfer Solution:** Optional standing belt attachments and caregiver foot positioning for added security.



Padded Lower Leg Supports



Ergonomic Base Design for Excellent Manoeuvrability

Specifications

	S332020
Overall Width	61cm
Overall Depth	64cm
Overall Height	118cm
Height of Base Plate	6cm
Base Plate Width (measured from rear)	51cm
Base Plate Width (measured from front)	36cm
Padded Leg Support (W x H)	16cm x 19cm
Padded Leg Support Height Adjustment	37cm - 54cm
Handle Height (from base plate - bottom/middle/top)	83cm/98cm/113cm
Unit Weight	18.5kg
Rear Castors Diameter	12cm
Safe Working Load (SWL)	150kg