

Presenting the Cobalt Air Chair Range

COBALT
HEALTH

Slimline Air Chair

Code: C01-SLIMLINE-BLU



Classic Air Chair

Code: C01-CLASSIC-BLU



Bariatric Air Chair

Code: C01-BARIATRIC-BLU



Our range of Air Chairs are specifically designed for use in all environments where users require high levels of comfort coupled with low – high risk pressure care.

Each Air Chair employs a combination of soft high-density foam, adjustable air cells in the seat and lumbar sections and a multi-directional stretch cover.

Every chair is highly adjustable, with independent tilt-in-space, recline and leg raise functions providing a wide array of repositioning options.

Additional features include drop-down armrests and fold-back head wings that can be lowered to allow for side transfers and to assist with cleaning.

A height-adjustable footrest, featuring independent flip up foot pads and an adjustable head pillow complete the comfort features of our versatile air chair range.

The Air Chair is available in Slimline, Classic and Bariatric, to cater for a wide range of users.

SWL Slimline & Classic	150kg
SWL Bariatric	200kg



**HELPING
CARE
HAPPEN**

Air Chair Features



Optional Tray Table



Footrests Included




Adjustable Air Cell



3 Brake & 1 Directional Castor

Slimline Air Chair



Safe Working Load	150kg
Overall Dimensions	130 x 65 x 125cm
Seat Width	45cm
Seat Depth	56cm
Seat Height	70cm
Seat Recline Angle	97-136°
Seat Tilt Angle	0-17°
Leg Rest Angle	0-90°
Order Code	C01-SLIMLINE-BLU

- Specialised pressure relief & support for smaller users
- Adjustable air cells in seat and lumbar sections
- Tilt-in-space, recline and leg raise functions
- 3 non-marking brake castors and 1 directional castor
- Fluid resistant, vapor permeable 4-way stretch cover
- Covers can be easily cleaned or replaced
- Fire retardant cover and foam
- Unit Weight 50kg
- 3-year warranty on frame, 12-month warranty on gas struts, upholstery, foam and cells

Classic Air Chair



Safe Working Load	150kg
Overall Dimensions	130 x 71 x 125cm
Seat Width	51cm
Seat Depth	56cm
Seat Height	70cm
Seat Recline Angle	97-136°
Seat Tilt Angle	0-17°
Leg Rest Angle	0-90°
Order Code	C01-CLASSIC-BLU

- Specialised pressure relief & support for most users
- Adjustable air cells in seat and lumbar sections
- Tilt-in-space, recline and leg raise functions
- 3 non-marking brake castors and 1 directional castor
- Fluid resistant, vapor permeable 4-way stretch cover
- Covers can be easily cleaned or replaced
- Fire retardant cover and foam
- Unit weight 50kg
- 3-year warranty on frame, 12-month warranty on gas struts, upholstery, foam and cells

Bariatric Air Chair



Safe Working Load	200kg
Overall Dimensions	130 x 80 x 125cm
Seat Width	60cm
Seat Depth	56cm
Seat Height	70cm
Seat Recline Angle	97-136°
Seat Tilt Angle	0-17°
Leg Rest Angle	0-90°
Order Code	C01-BARIATRIC-BLU

- Specialised pressure relief & support for bariatric users
- Adjustable air cells in seat and lumbar sections
- Tilt-in-space, recline and leg raise functions
- 3 non-marking brake castors and 1 directional castor
- Fluid resistant, vapor permeable 4-way stretch cover
- Covers can be easily cleaned or replaced
- Fire retardant cover and foam
- Unit weight 60kg
- 3-year warranty on frame, 12-month warranty on gas struts, upholstery, foam and cells

COBALT

HEALTH

FOR FURTHER INFORMATION OR DEMONSTRATIONS, PLEASE CONTACT US

Victoria

vicsales@cobalthhealth.com.au

(03) 8521 8444

43 Paramount Bvd, Cranbourne West VIC 3977

New South Wales

sales@cobalthhealth.com.au

(02) 8853 1100

34/9 Salisbury Rd, Castle Hill NSW 2154

cobalthhealth.com.au
1800 417 822