

Airline Travel and Battery Information

Power Chairs

Jazzy Carbon Ultra-lite

10.4AH (252Wh) Lithium-ion (standard)
20AH (504Wh) Lithium-ion (optional)

Airline Compliant
Not Compliant



Jazzy Passport

18.2AH, 24V Lithium-ion (436.8Wh) (standard)
10.4AH battery (optional)

Not Compliant
Airline Compliant



Scooters

Go Go Carbon

12AH, 24V Lithium-ion (288Wh)

Airline Compliant



Go Go Endurance Li

16AH Lithium-ion (2) x 8AH (385Wh) (standard)
(1) x 8AH (192Wh)
(Requires removal of x1 battery pack for compliance)

Not Compliant
Airline Compliant



Quest

(1) 11.5AH Lithium-ion (241.2Wh)

Airline Compliant



Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Pride lithium-ion battery are permitted.