

# EASY RIDER

## Scooter



“State of the art  
performance and design”



The Easy Rider has a striking new design,  
stunning instrumentation and offers  
precise handling that must be  
experienced to be appreciated.





# EASY RIDER

## Scooter

### Superb Comfort, Style & Performance...



Telescopic front motorcycle suspension for optimal handling



Motorcycle style wheels with chrome mud guards



Twin dial digital dash



All round lighting

- Top speed up to 10kph (6.2mph)\*.
- Maximum range up to 50km (31 miles) on full battery charge\*.
- Comfortable rotating sport style seat with sliding and reclining facility.
- Height adjustable headrest and angle adjustable armrests to suit user requirements.
- Angle adjustable handle bars allows comfortable positioning for the user (requires tools).
- Throttle lever with automatic safe-braking system.
- Additional Hand Brake for added security.
- Telescopic front motorbike suspension for optimal handling.
- Illuminated controls feature a twin dial digital dash.
- Trip computer with odometer for measuring distance travelled.
- Easy to use controls.
- All round lighting, including large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark.
- Twin adjustable rear view mirrors for extra visibility.
- Large accommodating footwell with anti-slip rubber mat.
- Motorbike style wheels with modern chrome mud guards.
- Rear bumper provides additional protection.
- Rear anti tip wheels as standard.
- Lap belt included for added safety.
- Large seat back storage pocket for keeping essentials safe when out and about.
- Anti-roll back system safety device prevents scooter from rolling backwards on hills.
- Freewheel facility allows movement of scooter without switching the motor on.
- Optional lockable and removable rear pannier box.
- Recommended MK Gel batteries: 2 x 12V (97Ah)\*\*.
- Weight capacity 180kg (28st/392lb).

\* Subject to conditions of use.  
\*\* Batteries not included

CODE: SR003BLKAU

DESCRIPTION: Easy Rider



Length:  
166cm (65")




Height:  
115cm (45")



Width:  
76cm (30")



Max Capacity:  
180kg (28st/392lb)



Ground Clearance:  
12cm (5")



Max Gradient:  
10°