



<15.4 km
Range per Charge



113 kg
Weight Capacity



61.5 cm
Turning Radius



Jazzy®
PASSPORT

Simply Fold and Travel

Get ready for the trip of a lifetime! The Jazzy® Passport brings you comfort, style, and convenience as you wait for your flight and explore your destination. This compact, lightweight power chair folds easily in just a few steps, making airport check-in and travel a breeze. With ample storage, the Jazzy® Passport has everything you need to embrace the adventures ahead. Experience the ease of lightweight travel on even terrain—wherever your journey takes you!



Jazzy[®] PASSPORT



Folds Easily for Effortless Transportation



Easy-Access Free-Wheel Levers



Under Seat Mesh Bag and Cup Holder

RANGE PER CHARGE (MAX) ¹	Up to 15.4 km
MAXIMUM WEIGHT CAPACITY	113 kg
TURNING RADIUS	61.5 cm
MAXIMUM SPEED ¹	5.7 km/h
LENGTH & WIDTH ²	L: 101.6 cm W: 59.6 cm
WEIGHT OF HEAVIEST PIECE ⁴	24.4 kg
GROUND CLEARANCE	7.6 cm at frame
SEAT SIZE WIDTH & DEPTH	W: 43 cm X D: 43 cm (17" X 17")

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Standard 18.2 Ah battery does not meet airline requirements. Optional 10.4 Ah battery makes this power chair airline compliant.
4. Seat weight included. Battery weight not included.

STANDARD ELECTRONICS	50A, PG VSI
DRIVE TRAIN	2 x 180W brushed motors, rear-wheel drive
BRAKING SYSTEM	Regenerative and electro-mechanical
SUSPENSION TYPE	Limited
FRONT CASTORS	8" solid
DRIVE WHEELS	12" pneumatic
BATTERY REQUIREMENTS ³	(1) 18.2Ah, 24V Lithium-ion (436.8Wh) (standard)
BATTERY CHARGER	Off-board, 2A Lithium-ion
BATTERY WEIGHT	3.9 kg
COLOUR OPTIONS	■ Pewter

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.



Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy Passport and its 10.4AH lithium-ion battery are permitted.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

20-24 Apollo Drive, Hallam, Victoria 3803 | (03) 8770 9600 | www.pridemobility.com.au



Quality
ISO 9001
SAI GLOBAL

AU_PASSPORT_16_01_25