

Air Comfort Compact

Comparison Guide

Air comfort

KCare
Healthcare Solutions



	Lift Chair Single Motor			Lift Chair Twin Motor			Mobile Lift Chair Single Motor			Mobile Lift Chair Twin Motor			Mobile Lift Chair Twin Motor - Titanium			
Size	Small	Medium	Large	Small	Medium	Large	Small	Medium	Large	Small	Medium	Large	Petite	Small	Medium	Large
Code	AC59018	AC59019	AC59020	AC59064	AC59065	AC59066	AC59098	AC59087	AC59099	AC59100	AC59090	AC59101	AC59124	AC59125	AC59126	AC59127
Material	VersoTela Fusion			VersoTela Fusion			VersoTela Fusion			VersoTela Fusion			VersoTela Theme			
Overall Width	750 mm	750 mm	820 mm	750 mm	750 mm	820 mm	750 mm	750 mm	820 mm	750 mm	750 mm	820 mm	750 mm	750 mm	750 mm	820 mm
Overall Height	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm	1060 mm
Seat Length	820 mm	820 mm	890 mm	820 mm	820 mm	890 mm	820 mm	820 mm	890 mm	820 mm	820 mm	890 mm	820 mm	820 mm	820 mm	890 mm
Seat Width	500 mm	500 mm	600 mm	500 mm	500 mm	600 mm	500 mm	500 mm	600 mm	500 mm	500 mm	600 mm	500 mm	500 mm	500 mm	600 mm
Seat Depth	440 mm	510 mm	510 mm	440 mm	510 mm	510 mm	440 mm	510 mm	510 mm	440 mm	510 mm	510 mm	410 mm	440 mm	510 mm	510 mm
Floor To Seat Height	480 mm	520 mm	520 mm	480 mm	520 mm	520 mm	480 mm	520 mm	520 mm	480 mm	520 mm	520 mm	480 mm	480 mm	520 mm	520 mm
Backrest Height	680 mm	740 mm	740 mm	680 mm	740 mm	740 mm	680 mm	740 mm	800 mm	680 mm	740 mm	800 mm	650 mm	680 mm	740 mm	800 mm
Max User Weight	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg	140 kg
Warranty	36 months - Frame/Foam 12 months - Gas Struts/Air Bags/Castors/Material			36 months - Frame/Foam 12 months - Gas Struts/Air Bags/Castors/Material			36 months - Frame/Foam 12 months - Gas Struts/Air Bags/Castors/Material			36 months - Frame/Foam 12 months - Gas Struts/Air Bags/Castors/Material			36 months - Frame/Foam 12 months - Gas Struts/Air Bags/Castors/Material			
Brand	Air Comfort			Air Comfort			Air Comfort			Air Comfort			Air Comfort			

Sales, administration and general enquiries National 1300 783 783 | customerservice@kcare.com.au | kcare.com.au

All information is considered to be true and correct at the date of publication. Changes in circumstances after the time of publication may impact on the accuracy of the information.

**Improving your
life everyday™**