

POWERSTROLL

Wheelchair Accessories

Power**STROLL**

S•Drive

*"Takes the strain
out of pushing
a wheelchair!"*



CE

This excellent Powerstroll will turn most manual wheelchairs into carer-controlled powered wheelchairs. It takes the strain out of pushing a heavy wheelchair and makes easy work of hills, ramps and uneven surfaces.

driveTM
www.drive-medical.com.au

The S-Drive Powerstroll can easily be fitted and removed within seconds through its unique one-click docking station (patented) and the new push-and-slide battery removal system.



Easy to fit & remove



Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control



Fits most standard wheelchairs up to 50cm (20") seat width

Superb Comfort, Style & Performance...

- » Allows speed of up to 6.5kph (4mph)*.
- » Travels up to 16km (10 miles) on a full battery charge*.
- » The powerstroll can easily be fitted or removed within seconds.
- » The excellent PG S-Drives controller features a fully programmable and comprehensive fault management system.
- » New push-and-slide battery removal system.
- » It is very portable, making transportation simple.
- » Fits most standard wheelchairs up to 50cm (20") seat width.
- » Additional parts available for wider wheelchairs.
- » Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control.
- » Unique one-click docking station for ease of assembly (patented).
- » Comes complete with a carry holdall with shoulder strap.
- » Simple to manoeuvre and reverse.
- » The powerstroll has a powerful 150 watt (200 watt-HD), 24V motor, capable of coping with most terrains.
- » Maintenance free, airline friendly 24v 12Ah valve regulated sealed lead acid battery pack.
- » 2A charger as standard for off-board charging.
- » Weight capacity 135kg (21st/300lb), this includes the wheelchair and user combined.
- » Heavy Duty model - maximum weight capacity 180kg (28st/400lb).

* Subject to conditions of use.

CODE: **PWCPP009** (Standard)

CODE: **PWCPP009HD** (Heavy Duty)

DESCRIPTION: **S-Drive Powerstroll**



Total Weight:
19.2kg (42lb)



Weight Without Battery:
10.2kg (22.5lb)



Weight Battery Pack Only:
9kg (19.5lb)



Max Range:
16km (10 Miles)*



Max Gradient:
8°



Max Weight Capacity:
135kg (21st/300lb) - ST
180kg (28st/400lb) - HD

Leaflet Code: **LLO01**

If you would like to view this product please contact your local stockist, details below.