

END-TO-END | Scooter Specifications

																						
	QUEST®	GO GO® CARBON	GO GO® e	GO GO® ULTRA X	GO GO® ELITE TRAVELLER 2	GO GO® LX W/ CTS SUSPENSION	GO GO ENDURANCE Li®	ZERO TURN 8	BAJA™ BANDIT	ZERO TURN 10	PATHRIDER® ES10	PATHRIDER® TWIN	PATHRIDER® 10 2.0	PATHRIDER® 10 DX	VICTORY® 10 LX W/ CTS SUSPENSION	MAXIMA	PATHRIDER® ENDURANCE	CELEBRITY MFP	PATHRIDER® 130 XL	PATHRIDER® 140 XL	OUTBACK®	SPORT RIDER
RANGE PER CHARGE (MAX) ¹	Up to 19.3 km	Up to 16 km	Up to 14.4 km	Up to 11 km	Up to 16.7 km	Up to 16.9 km	Up to 23 km	Up to 16 km	Up to 23.3 km	Up to 28.8 km	Up to 17.9 km	Up to 25 km	Up to 18 km	Up to 25 km	Up to 26 km	Up to 31.4 km	Up to 54.7 km	Up to 30 km	Up to 25 km	Up to 35 km	Up to 35 km	Up to 35 km
MAXIMUM WEIGHT CAPACITY	113 kg	136 kg	120 kg	118 kg	136 kg	136 kg	147 kg	136 kg	181 kg	170 kg	158 kg	159 kg	181 kg	158.7 kg	181 kg	227 kg	158 kg	226 kg	158.7 kg	181 kg	158 kg	181 kg
TURNING RADIUS	82 cm	134.6 cm	125 cm	100.3 cm	94 cm	116.2 cm	122.5 cm	97.2 cm	134.7 cm	109 cm	132 cm	107 cm	132.5 cm	133.3 cm	152 cm	3W: 124.5 cm 4W: 134.6 cm	176 cm	158 cm	170 cm	210 cm right 208 cm left	243 cm	158 cm
MAXIMUM SPEED ¹	Up to 6 km/h	Up to 6 km/h	Up to 6.4 km/h	Up to 6.4 km/h	Up to 6.8 km/h	Up to 7 km/h	Up to 6.7 km/h	Up to 6.4 km/h	Up to 10 km/h	Up to 10 km/h	Up to 8.8 km/h	6.4 km/h	Up to 8.4km/h	Up to 8.9 km/h	Up to 9 km/h	Up to 9.3 km/h	Up to 9.7 km/h	Up to 10 km/h	Up to 10 km/h	Up to 10 km/h	Up to 10 km/h	Up to 10 km/h
LENGTH & WIDTH	L: 92.7 cm W: 49.5 cm ²	L: 102.3 cm W: 44 cm	L: 100 cm W: 49 cm ²	L: 100.3 cm W: 49.5 cm ²	L: 103.6 cm W: 49.4 cm ²	L: 104 cm W: 52.1 cm ²	L: 107 cm W: 52 cm ²	L: 116 cm W: 50.8 cm ²	L: 126 cm W: 56 cm ²	L: 122 cm W: 61.6 cm ²	L: 115.6 cm W: 57.1 cm ²	L: 116 cm W: 56 cm ²	L:120 cm W: 56.5 cm ²	L: 119.4 cm W: 56.2 cm ²	L: 119.4 cm W: 56.5 cm ²	L: 120.7 cm W: 68.6 cm ²	L: 146 cm W: 75.3 cm ²	L: 145 cm W: 70.7 cm ²	L: 136 cm W: 63.6 cm ²	L: 142.8 cm W: 69.2 cm ²	L: 156 cm W: 83 cm ²	L: 165 cm W: 78 cm ²
SEAT TYPE WIDTH & DEPTH (STANDARD)	Compact, foldable, molded plastic seat W: 41 cm (16") D: 33 cm (13")	Compact, foldable seat W: 35.6 cm (14") D: 27.5 cm (11")	Molded plastic, black vinyl W: 41.9 cm (16.5") D: 33 cm (13")	Compact, foldable seat W: 43 cm (17") D: 43 cm (17")	Compact, foldable seat W: 43 cm (17") D: 43 cm (17")	Compact, foldable seat W: 43 cm (17") D: 43 cm (17")	Compact, foldable seat W: 43 cm (17") D: 43 cm (17") or W: 50.8 cm (20") D: 45.7 cm (18")	Compact, foldable seat W: 43 cm (17") D: 43 cm (17")	Compact, foldable seat W: 46.5 cm (18.3") D: 44 cm (17.3")	Deluxe, pillow-top, high-back seat W: 49.5 cm (19.5") D: 44.5 cm (17.5")	Compact, foldable seat W: 46 cm (18") D: 43 cm (17")	High-back seat W: 48.3 cm (19") D: 55.2 cm 21.75")	High-back seat W: 43.2 cm (17") D: 43.2 cm (17")	Pillow-top, high-back seat W: 44.2 cm (17.5") D: 43 cm (17")	Pillow-top, high-back seat W: 46 cm (18") D: 46 cm (18")	High-back seat W: 55.9 cm (22") D: 50.8 cm (20")	High-back seat W: 55.9 cm (53") D: 50.8 cm (45.7")	Deluxe, pillow-top, high-back seat W: 49.5 cm (19.5") D: 44.5 cm (17.5")	Deluxe, pillow-top, high-back seat W: 49.5 cm (19.5") D: 44.5 cm (17.5")	Deluxe, pillow-top, high-back seat W: 49.5 cm (19.5") D: 44.5 cm (17.5")	Deluxe, pillow-top, high-back seat W: 49.5 cm (19.5") D: 44.5 cm (17.5")	Captains high-back, adjustable seat W: 51 cm (20") D: 56 cm (22")
SEAT-TO-DECK HEIGHT RANGE	31 cm fixed	43 cm fixed	41.9 cm - 45.7 cm	40 cm - 43.8 cm	41.6 cm - 45.3 cm	41.3 cm - 43.8 cm	41.2 cm x 45 cm	40.6 cm - 45.6 cm	47.5 cm suspension seat post	47.5 cm - 52.5 cm	37 cm - 42 cm	48.3 cm - 53.3 cm	43.8 cm - 47.6 cm	43.8 cm - 47.6 cm	46 cm - 49 cm	43.5 cm - 48.6 cm	45.5 cm - 53 cm	45.2 cm - 49 cm	47 cm suspension seat post	47 cm fixed	50 cm - 53 cm	48 cm fixed
SEAT-TO-GROUND HEIGHT RANGE	51 cm fixed	57 cm fixed	54.6 cm - 59.7 cm	55 cm - 60 cm	46.3 cm - 51.4 cm	52.7 cm - 56.5 cm	53.9 cm x 59 cm	53.3 cm - 60 cm	66 cm suspension seat post	58.4 cm - 63.5 cm	56.7 cm suspension seat post	63 cm - 67 cm	60 cm - 63.8 cm	60 cm - 63.8 cm	60 cm - 64.8 cm	56.1 cm - 61.2 cm	55.6 cm - 63 cm	66.7 cm - 70.5 cm	56.7 cm suspension seat post	66.7 cm - 70.5 cm	76 cm - 80 cm	24.8 cm fixed
GROUND CLEARANCE	5.8 cm at frame ³	3.8 cm at motor ³	5 cm at transaxle ³	6.3 cm at midframe 3.8 cm at motor ³	3.4 cm at transaxle ³	6.4 cm at frame 3.8 cm at motor ³	3.6 cm at motor ³	8.9 cm at frame 3.3 cm at motor ³	4 cm at motor skid plate	5 cm at motor ³	8.8 cm at frame ³	5.7 cm at motor ³	6 cm at motor ³	8.25 cm at frame ³	5.6 cm at motor ³	5.7 cm at motor ³	9.6 cm at transaxle ³	8.8 cm at motor	9 cm at transaxle ³	13.3 cm at midframe	14 cm at midframe	12.7 cm at shroud ³
OVERALL WEIGHT (NOT INC. BATTERY)	23.5 kg ⁴	16 kg	32.5 kg ⁴	36.3 kg ⁴	37.6 kg ⁴	39.4 kg ⁴	42.4 kg ⁴	53.9 kg ⁴	72.3 kg ⁴	82 kg ⁴	59.4 kg ⁴	68 kg ⁴	63 kg ⁴	84.8 kg ⁴	81.2 kg ⁴	3W: 83.2 kg 4W: 89.3 kg ⁴	100 kg ⁴	104.3 kg ⁴	99.8 kg ⁴	117.7 kg ⁴	124 kg ⁴	115 kg ⁴
WEIGHT OF HEAVIEST PIECE	23.5 kg (does not disassemble)	16 kg (does not disassemble)	31 kg (does not disassemble)	15.4 kg (front sections)	15.8 kg (front section)	16.6 kg (front section)	18.5 kg (front section)	21.8 kg (rear section)	31.8 kg (front section)	34 kg (rear section)	25.8 kg (front section)	27.3 kg (front section)	29.5 kg (front section)	29.5 kg (front section)	35.2 kg (rear section)	3W: 61.7 kg 4W: 67.3 kg (main frame does not disassemble)	34 kg (rear section)	94 kg (main frame does not disassemble)	34 kg (rear section)	93 kg (main frame does not disassemble)	99 kg (main frame does not disassemble)	92.3 kg (main frame does not disassemble)
SUSPENSION	N/A	N/A	N/A	N/A	N/A	CTS Front & Rear	CTS Front & Rear	CTS Front & Rear	Front, Rear & Seat Post	Front & Rear	N/A	Rear	N/A	Front & Rear	Front & Rear	N/A	Front & Rear	Front & Rear	Front & Rear	Front, Rear & Seat Post	Front & Rear	Front & Rear
3W/4W OPTION	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	4W	3W & 4W	4W	4W	4W	4W	4W	3W
FRONT TYRES	7" solid	7" solid	8" solid	7" solid	7" solid	7" solid	8" solid	7" solid	9" pneumatic	9" solid	10.4" solid	9" solid	10" pneumatic	10.4" solid	10" pneumatic	10" solid	13" pneumatic	13" pneumatic	13" pneumatic	13" pneumatic	14.5" pneumatic	19" pneumatic
REAR TYRES	7.8" solid	8" solid	8" solid	8" solid	8" solid	8" solid	8" solid	9" solid	9" pneumatic	10.75" solid	10.75" solid	10.75" solid	10" pneumatic	10.75" solid	10" pneumatic	10.75" solid	13" pneumatic	13" pneumatic	13" pneumatic	13" pneumatic	14.5" pneumatic	19" pneumatic
BATTERY REQUIREMENTS	(1) 11.5AH (lithium-ion, 241.2Wh) ⁹	(1) 12AH, 24V (lithium-ion, 288Wh) ⁹	(2) 12AH ⁵	(2) 12AH ⁵	(2) 18AH ⁵	(2) 18AH ⁵	16AH (lithium-ion, 2 x 8AH, 385Wh) ⁷	(2) 22AH ⁵	(1) 20AH, 36V (lithium-ion, 768Wh) ⁸	(2) 40AH ⁵	(2) 31AH ⁵	(2) 35AH ⁵	(2) 31AH ⁵	(2) 31AH (standard) or (2) 40AH (optional extra) ⁵	(2) 40AH ⁵	(2) 50AH ⁵	(2) 50AH (lithium-ion, 1280Wh) ⁸	(2) 70AH ⁵	(2) 55AH ⁵	(2) 70AH ⁵	(2) 100AH ⁵	(2) 70AH ⁵
BATTERY CHARGER	Off-Board 2.5A (lithium-ion)	Off-board, 2A (lithium-ion)	Off-Board, 2A	Off-Board, 2A	Off-Board, 2A	Off-Board, 2A	Off-board, 3.5A (lithium-ion)	Off-Board, 3.5A	Off-Board, 4A (lithium-ion)	Off-Board, 3.5A	Off-Board, 3A	Off-board, 3.5A	Off-Board, 3.5A (32AH)	Off-Board, 3.5A (32AH) or Off-Board, 5A (40AH)	Off-Board, 5A	Off-board, 5A	Off-Board, 8A (lithium-ion)	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A	Off-Board, 8A
BATTERY WEIGHT	2.3 kg	1.8 kg	2.5 kg each	4 kg each	6 kg each	6 kg each	5.4 kg each	6 kg each	9.5 kg	14 kg each	10.5 kg each	10.5 kg each	10.3 kg each	10.5 kg (standard) or 14 kg	14 kg each	16.6 kg	10.8 kg each	22 kg each	16.8 kg each	23.5 kg each	27.9 kg each	23.7 kg
COLOUR OPTIONS	White	Black, Yellow	Matte Black	Red	Red, Black, Blue, White ⁶	Red, Blue ⁶	Garnet Red, Euro Grey ⁶	Fire Opal, Pearl White, Sapphire Blue	Sage, Tan	Candy Apple Red	Red	Red, Blue	Red	Viper Blue	Red	Garnet Red, Ocean Blue	Black, Orange	Red, Blue	Candy Apple Red, Mercury Silver	Candy Apple Red, Yellow, Silver	Blue	Black

1 Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition **2** Width of base (not including seat) **3** Does not account for anti-tip clearance **4** Includes standard seat weight **5** Deep cycle (AGM or gel cell type recommended) **6** All shroud panels included **7** Standard 16AH battery does not meet airline requirements To make this scooter airline compliant please remove one of the 8AH batteries before travelling on an aircraft **8** Not suitable for air travel **9** Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy Carbon® Ultra-lite and its lithium-ion battery are permitted. **NOTE:** All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%

