

2700 Six Channel Cordless Falls Monitor

cura¹
Helping To Care



... cordless without compromise

2700 Six Channel Cordless Falls Monitor

No more cords between monitor and pads! Our modern, ergonomic cordless systems support a professional, clutter-free environment and fast, economical setup. This full-featured cordless monitor is our flagship monitor.

- Can be used with up to six cordless pads, mats or other triggering devices and one corded (premium) device.
- Acts as a room controller within a bed management system.
- Fault supervision for pressure pads (low battery, device lost).
- Battery operated and / or optional external plug pack.
- Comes with batteries, soft monitor cover and reset key.
- Up to 6 months battery life (nominal).
- All safety alerts sound even if alarm volume set to minimum, including low battery, AC lost, pad signal lost & nurse call lost.



WHAT'S IN THE BOX?



2700 Six Channel Cordless Falls Monitor

2742 Soft Monitor Cover

Batteries

2642 Caregiver Reset Key

QuickStart User Guide

COMPATIBLE PRODUCTS

- All cura¹ cordless pads, mats, sensors & buttons.
- All cura¹ premium pads & mats with RJ10 plugs.
- All cura¹ universal pads & mats.
- 2406 Plug Pack 6 V.



Certifications may vary between products.

cura¹
Helping To Care

10 Gray Street, Kilkenny, South Australia 5009

sales@cura1.com | www.cura1.com

+61 (0) 8 8268 8343 (International)

1300 12-CARE (2273)