

CONFIGURA BARIATRIC ASSESSMENT FORM

Prescriber's Name: _____ Client's Name: _____
 Number: _____ Address: _____
 Email: _____ Number: _____
 Funding Body: _____



Size Areas highlighted in grey are standard specification

Tilt-In-Space Configura Bariatric (Dual Motor with Independent Legrest) in Black

Product Code	Seat Height	Seat Width	Seat Depth	SWL
CR5470	406 mm (16") <input checked="" type="checkbox"/>	609 mm (24") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	254kg
CR5471		660 mm (26") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	
CR5472		711 mm (28") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	
CR5473	457 mm (18") <input type="checkbox"/>	609 mm (24") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	254kg
CR5474		660 mm (26") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	
CR5475		711 mm (28") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	
CR5476	508 mm (20") <input type="checkbox"/>	609 mm (24") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	317kg
CR5477		660 mm (26") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	
CR5478		711 mm (28") <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	

Standard Configura Bariatric (Dual Motor) in Duratek

Product Code	Seat Height	Seat Width	Seat Depth	SWL
CR5412	406 mm (16") <input type="checkbox"/>	609 mm (24") # ² <input type="checkbox"/> 660 mm (26") <input type="checkbox"/> 711 mm (28") # ³ <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	254kg
CR5414		609 mm (24") # ² <input type="checkbox"/> 660 mm (26") <input type="checkbox"/> 711 mm (28") # ³ <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	
CR5420		762 mm (30") # ² <input type="checkbox"/> 812 mm (32") <input type="checkbox"/> 863 mm (34") # ³ <input type="checkbox"/>	457 mm (18") # ¹ <input type="checkbox"/> 508 mm (20") <input type="checkbox"/> 558 mm (22") <input type="checkbox"/>	

(#1) 50mm (2") Depth Adjustment Kit <input type="checkbox"/>	(#2) 24"/ 30" Seat Width <input type="checkbox"/>	(#3) 28"/ 34" Seat Width <input type="checkbox"/>
<input type="checkbox"/> Duratek (Standard) CA2432	<input type="checkbox"/> Black (TIS) CA2801	Configura Bariatric Armrest, 24" (Std TIS), 30" (Supa), Duratek, PAIR CA2514 Swap the 26"/32" LHS & RHS armrests around to create a seat width of 28"/34". The overhang will now be on the outside of the chair.

Extra	Posture
Channel Legrest <input type="checkbox"/>	
<input type="checkbox"/> Duratek CA2446 <input type="checkbox"/> Charcoal VP CA2544 <input type="checkbox"/> Black VP CA2802	<input type="checkbox"/> Waterfall Backrest (std) 
<input type="checkbox"/> Configura Battery Backup CA2402	<input type="checkbox"/> Laterals CA2413  SOLD IN A PAIR
	<input type="checkbox"/> Profiled Headrest  <input type="checkbox"/> Charcoal VP (Standard) CA2448 <input type="checkbox"/> Black VP (TIS) CA2800
Pressure	
Rise Angle Setting	
20 degrees of tilt (std) <input type="checkbox"/>	<input type="checkbox"/> Visco a Cushion (std) 
30 degrees of tilt (TIS model only) <input type="checkbox"/>	<input type="checkbox"/> Castellated Foam 
Handset Control Position	
RHS (std) <input type="checkbox"/>	<input type="checkbox"/> Charcoal VP Cover (Standard) CA2434
LHS <input type="checkbox"/>	<input type="checkbox"/> Black VP Cover (TIS) CA2804
	<input type="checkbox"/> Cushionair Alternating Seat Cushion w/ Pump 
	<input type="checkbox"/> Charcoal VP Cover (Standard) CA2720 <input type="checkbox"/> Black VP Cover (TIS) CA2725

Chair Assessment and Delivery Environmental Questionnaire

When carrying out an assessment please fill out the below form prior to ordering the chair.

Please give as much detail as possible.

ACCESS

Is there parking/access for a large van? _____

Are there steps to the property? _____

Is it easy access to the property? _____

INTERNAL

Are the doors and corridors wide enough for a chair to go through? Min. width for most chairs is 650mm.

Does the chair have to go upstairs? If yes, please specify number of steps and if a stair lift is fitted.

Is space being cleared for the chair? _____

Is there a suitable power supply for the chair? Yes No

INFECTION CONTROL

Are there pets? _____ Is the house smoke free? _____

OTHER

Does the driver need to be aware of any other points? _____