

Corzino Stroller

The Patron Corzino Stroller is a lightweight, easy-folding adaptive stroller designed for short-term use, helping to make everyday outings easier. Whether it's a quick trip to the shops, a family adventure or navigating busy streets, the Corzino provides a safe and comfortable ride for children who need occasional mobility support. Available in four sizes, it's a compact, stylish and practical solution for families looking for a reliable stroller that fits easily into daily life.

Compact, adaptable, and built for everyday adventures

Ideal for short-term use and fatigue management, enabling safe community access

Four size options (30, 34, 38, 42 cm seat widths) to suit children through to teens

Lightweight alloy frame for effortless manoeuvrability



High weight capacity (up to 75 kg for sizes 38, 42) with an exceptional size-to-weight ratio

Anti-escape options to support children with behaviours of concern in public settings

Great for urban conditions, public transport, and travel



5-point comfort harness for enhanced safety and support



Adjustable seat depth, footrest height, and backrest height extension (size 38, 42 only) for long-term adaptability



Ultra-compact umbrella folding design for easy storage and transport



Large solid wheels (175/265 mm) for smooth navigation over urban and light off-road terrain

ACCESSORIES

Extended Side Protection

Small (Corzino 30):
2510-6080-000
Medium (Corzino 34):
2510-6072-000
Large (Corzino 38/42):
2510-6073-000



Sun Canopy

2510-7519-000



Rain Cover

(Must be used with Sun Canopy).
For Corzino 30/34:
2510-7522-000
For Corzino 38/42:
2510-7523-000



Upper Limb Support

2510-6140-000



Storage Basket

Small (Corzino 30/34):
2510-8010-100
Large (Corzino 38/42):
2510-8010-200

Pen Key Release Hip Belt

X-Small (30cm Pad):
1405-6210-025
Small (36cm Pad):
1405-6211-025
Medium (46cm Pad):
1405-6212-025



Package Includes:

- Adjustable seat depth
- Soft and padded upholstery
- Removable upholstery
- Comfortable 5-point support harness
- Comfort side protection
- Back support recline adjustment (2 positions)
- Height-adjustable and foldable foot supports
- Solid tyres
- Removable wheels
- Heel strap

(Storage Basket not included).



Size	Code	Approx. Age of User	Weight Limit	Seat Width	Seat Depth	Back Support Height	Back Support Angle	Footplate Height	Approx. Stroller Weight
30	2510-0001-030	4-8 years	45kg	30cm	28-34cm	55cm	Fixed	27-41cm	13kg
34	2510-0001-034	6-10 years	55kg	34cm	29-35cm	67cm	Fixed	29-43cm	13.1kg
38	2510-0001-038	8-12 years	75kg	38cm	32-39cm	65cm (85cm with extension - supplied as Std.)	Fixed	30-47cm	15kg
42	2510-0001-042	10-14 years	75kg	42cm	33-41cm	69cm (88cm with extension - supplied as Std.)	Fixed	32-48cm	16kg