

Hero-XL Car Seat

The introduction of the Hero-XL – Australasia's first extended rear-facing car seat for users with disabilities – means children can remain rear-facing for longer while getting the postural support they need.

Super Support!

Safety First

From the integrated 5-point harness with quick release, to the dynamic headrest and side impact protection, everything about the Hero-XL was designed with safety in mind.

Fully Adjustable

Tailor a unique fit for your unique child, with adjustable lateral trunk supports, tray tables and medial thigh supports, and customisable seat depth.

Years of Support


Your Hero-XL rear-facing car seat will serve your child for years, with the ability to accommodate users up to 125cm in height and 36kg in weight.

Easy Installation

The Hero-XL does not require an ISOFIX connection and can be quickly secured via your vehicle's seatbelt system. Includes straps to lower tether points.



Rear-facing seats can protect children from the forces of braking and deceleration, allowing those with limited head and trunk control to travel safely in comfort.



The Hero-XL includes a range of features to enhance that comfort further, with soft, padded tray tables, adjustable seat depth, removable medial thigh support (to prevent hip abduction) and lateral trunk supports which can be adjusted for weight, height and angle.

ACCESSORIES

Seat Depth Adaptor
4182-0280-000



Medial Thigh Support
4182-1510-000



Combination Medial Thigh Support and Depth Adaptor
4182-0895-000



Tray Table
4182-5544-000



Additional Lower Tether Anchorage Straps
4182-6953-000



Seat Protector Mat
4182-0190-000
(to go underneath car seat)

Package includes

Hero-XL extended rear facing special purpose car seat in Grey/Black complete with:

- Accessory mounting frame
- Adjustable lateral supports
- Installation straps for lower tether points

Code	Seat Width (knee)	Seat Width (pelvis)	Seat Depth	Back Height	Distance Between Laterals	Total Width	Capacity
4182-0010-090	27cm	24cm	25-34cm	28-45cm	13-20cm	43cm	36kg or 125cm Tall