

Piper Comfort Stroller

The Piper Comfort Stroller is a lightweight, low-support adaptive stroller that adapts to the needs of growing children. Compact and foldable, it's perfect for families on the go, effortlessly navigating tight spaces and public transport while providing a comfortable and secure ride. With its sturdy construction and compatibility with various accessories, the Piper Comfort is designed for everyday adventures, offering the support children need to thrive in their community.

Reliable support for growing explorers

Low seat height for seamless transfers

Designed to accommodate growing children for prolonged use

5-point harness for improved safety and support

Anti-escape options to support children with behaviours of concern in public settings

Durable, low-maintenance tyres for hassle-free use

Carry basket included for convenience



Single-point release harness for quick and easy removal



Adjustable and fully enclosed back support recline for enhanced comfort and security



Ultra-compact folding design for effortless storage and transport



Large solid wheels (175/265 mm) for smooth rides across urban & light outdoor terrain

ACCESSORIES

Seat Inlay
2510-6080-000



Upper Limb Support
2510-6140-000




Footplate Height Reducing Spacer
(Reduces footplate height by 5 or 10cm).
2510-6170-000



Sun Canopy
2510-7519-000



Rain Cover
(Must be used with Sun Canopy).
2510-7522-000



Pen Key Release Hip Belt
(30cm Pad).
1405-6210-025



Package Includes:

- Removable upholstery
- Comfortable 5-point support harness
- Back support recline adjustment
- Height-adjustable and detachable 1 piece foot support
- Solid tyres
- Removable wheels
- Swivel locks on front wheels
- Basket under stroller for storage



Code	Weight Limit	Seat Width	Seat Depth	Back Support Height	Back Support Angle	Footplate Height	Approx. Stroller Weight
2510-0003-034	50kg	23/27/34cm*	24/29cm*	68cm	90-135°	10 – 30	14.5kg

*full range requires use of Seat Inlay accessory