

Ella Stander

From early paediatrics to adults, the Ella Stander adapts and grows with the user, delivering tailored support and positioning in either prone or supine.

The stander that does it all

Stepless Support Adjustments

Fine-tune pelvic, trunk, and head supports to meet complex postural needs.

Removable, Wipeable Upholstery

Easy-care materials designed for real-world use and hygiene needs.

Future-Proof with Growth Accommodation

Designed to support users across many years with a huge growth range.



Folding Base for Compact Storage

Disassembles easily with just two screws; ideal for transport between settings.

Advanced Foot Sandal Design

Secure, fully angle-adjustable foot supports with high edges for heel stability and no slip-offs.

Independent Knee and Foot Adjustments

Adapt angles in all directions, including abduction, height, width, and depth, for tailored leg and foot positioning.



All-in-One Supine and Prone Standing

Switch between prone and supine standing in under two minutes using a single tool.



Prone set-up

Supine set-up

Two Choices of Base

Offering 90° Tilt Facilitate safe, effortless horizontal-to-vertical transfers with either powered or manual control.



Four Package Options to Suit Clinical Need

Choose between Standard or Complex setups, with essential accessories included, plus a choice of manual or powered base. Alternatively you can 'Build Your Own' package.

Standard Package includes:

- Mobile base with powered or manual tilt
- Large capacity abduction kit
- Lockable castors
- Adjustable knee and foot supports
- Pelvic and trunk support
- Belt
- Tray
- Swing-away straight laterals
- Headrest
- Tilt angle indicator and measurement guide
- Navy wipe-clean upholstery

Complex Package Includes:

- All Standard features, plus:
- De-rotational pelvic and trunk supports
 - Flexi-laterals for trunk
 - Shoulder protraction pads
 - Harness fitting kit
 - Rotational adjustment knee support

Size	Code	Max User Weight	Approx User Height	Footplate to Axilla	Chest / Hip Width	Foot to Patella	Foot to Pelvis	Abduction (Max Angle)	Kneepads to Footrest
1	Standard: 5221-0020-000	30kg	60-95cm	43-68cm	12-20cm	14-45cm	30-54cm	67°	0-25°
2	Standard: 5222-0020-000 Complex: 5222-0030-000	40kg	90-130cm	57-97.5cm	17-25cm	18-55cm	47.5-78.5cm	45°	0-25°
3	Standard: 5223-0020-000 Complex: 5223-0030-000	55kg	120-155cm	80-128cm	19-30cm	20-65cm	50-104cm	43°	0-25°
4	Standard: 5224-0020-000 Complex: 5224-0030-000	80kg	130-180cm	90-150cm	25-41cm	25-85cm	62.5-122.5cm	55°: 130cm* 35°: 180cm*	0-25°

*User Height

Accessories/Build Your Own

Base

Size 1, Manual 5220-0111-000
 Size 2, Manual 5220-0112-000
 Size 3, Manual 5220-0113-000
 Size 1, Power 5220-0121-000
 Size 2, Power 5220-0122-000
 Size 3, Power 5220-0123-000
 Size 4, Power 5220-0124-000

Foot Supports (Pair)

Size 1 5220-6071-000
 Size 2 5220-6072-000
 Size 3 5220-6073-000



Knee Supports (Pair)

Size 1 5220-1761-000
 Size 2 5220-1762-000
 Size 3 5220-1763-000



Knee Support Rotation - Module

5220-7200-000



Pelvic Support Boards

Size 1, Standard 5220-1911-000
 Size 2, Standard 5220-1912-000
 Size 3, Standard 5220-1913-000
 Size 4, Standard 5220-1914-000
 Size 1, Rotational 5220-1941-000
 Size 2, Rotational 5220-1942-000
 Size 3, Rotational 5220-1943-000
 Size 4, Rotational 5220-1944-000

Trunk Support Boards

Size 1, Standard 5220-2911-000
 Size 2, Standard 5220-2912-000
 Size 3, Standard 5220-2913-000
 Size 4, Standard 5220-2914-000
 Size 1, Rotational 5220-2961-000
 Size 2, Rotational 5220-2962-000
 Size 3, Rotational 5220-2963-000
 Size 4, Rotational 5220-2964-000

Rotational Lumbar Support Pad

S 5220-2601-000
 M 5220-2602-000
 L 5220-2603-000



Support Belts

Size 1 5220-6821-000
 Size 2 5220-6822-000
 Size 3 5220-6823-000
 Size 4 5220-6824-000



Swing-Away Lateral Supports (Pair)

Size 2 5220-3321-000
 Size 3 & 4 5220-3331-000
 Size 3 & 4 Inset 5220-3333-000



Flexi Pad, Swing-Away Lateral Supports (Pair)

Size 2 5220-3421-000
 Size 3 & 4 5220-3431-000
 Size 3 & 4 Inset 5220-3430-000



Elbow Blocks

Small 5220-7351-000
 Large 5220-7352-000

Standard Tray

For Size 1 5220-7121-000
 For Size 2 5220-7122-000
 For Size 3 & 4 5220-7123-000



Tray with Bowl

For Size 1 5220-7181-000
 For Size 2 5220-7182-000
 For Size 3 & 4 5220-7183-000



Centrepont Harness

X-Small 1409-6660-017
 Small 1409-6661-017
 Medium 1409-6662-017
 Large 1409-6663-017



Shoulder Protraction Pads (pair)

Size 1 5220-2751-000
 Size 2 5220-2752-000



Harness Fitting Kit

5220-6150-000



Adjust-A-Plush Head Support

Standard 1480-4330-281
 Medium 1480-4320-281
 Narrow 1480-4310-281



Plush Head Support

10" 5220-4610-000
 14" 5220-4614-000

Tool Free Kit

5220-8160-000