

# ROYALE

## ROBOOTER

Smart style. Everyday simplicity.

● E40



The Royale Robooter E40 Powerchair combines sleek design with practical performance for everyday use. Lightweight, foldable and stylish, it's ideal for home and community mobility where comfort and control matter most.

# E40

## Sleek powerchair with folding design

Modern comfort meets everyday convenience – for those who want powered mobility that suits their life and routine, seamlessly.

Compact yet powerful, the E40 is built for effortless mobility. A slim white frame and smooth joystick control make it easy to navigate hallways, cafés or apartments, while the quick-fold design allows for easy packing in the car boot within seconds.

With attendant controlled remote included, the E40 makes powered mobility simple, portable and easy to live with.



### Effortless folding frame

Stores easily with no disassembly



### Lightweight 23kg design

Ultra compact for transport



### Compact turning circle

Small footprint with exceptional manoeuvrability







## E40

- ▶ **Flip up armrests**  
Safe transfers
- ▶ **Ergonomic seat & backrest**  
Comfort for daily use
- ▶ **LED front & rear lights**  
Improved visibility

### Attendant control

Via app or bluetooth remote




### Storage included

Rear backpack & underseat basket



### Travel up to 23km

With 20 AH battery



## Specifications

Frame Colour	White
Frame Material	Aluminium
Safe Working Load (SWL)	150 kg
Battery Type	Lithium-Ion
Battery Size	Standard: 20 AH Optional: 12 AH or 25 AH
Range	17 km (12 AH) 23 km (20 AH) 29 km (25 AH)
Max Speed	8 km/h
Size (W x L x H Assembled)	62.4 x 100 x 93 cm
Seat Width	42 cm

Seat Depth	43 cm
Backrest Height	46 cm
Backrest Angle	18°
Distance Between Armrests	48 cm
Weight (Without Battery)	23 kg
Front Wheels	8" PU
Rear Wheels	10" pneumatic
Max Safe Incline Slope	6°
Turning Circle	90 cm
Clearance	4 cm
Folded Size (W x L x H)	62.4 x 70 x 45 cm