

EDGE 3

EDGE 3
shown in Green Machine with
iLevel® power seat elevation and
Q-Logic 3 with Bluetooth®



MADE FOR EVERYDAY PERFORMANCE

The Edge® 3 is built to deliver a smooth, stable ride and responsive performance for everyday use. Powerful motors and Smooth Ride Suspension (SRS) provide exceptional comfort and handling across a range of environments. Optional iLevel® technology offers up to 30.5 cm of power adjustable seat elevation at walking speed (7.2 km/h), supporting greater independence throughout the day. Standard features include front and rear LED fender lights, a USB charger, and bold drive wheel accents.



QUANTUM®



quantumrehab.com.au

Specifications

DRIVE WHEELS (Mid-Wheel Drive)	35.5cm (14") Pneumatic or 14x4" wide Flat Free Tyres
CASTER WHEELS	Front: 15.2 cm (6") solid Rear: 15.2 cm (6") solid
ANTI -TRIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 10 km/h Up to 7.2 km/h w/ iLevel activated
GROUND CLEARANCE	6.8 cm (2.7")(battery Tray)
TURNING RADIUS	52 cm (20.5")
OVERALL LENGTH ²	90 cm (35.5") without front riggings
BASE WIDTH	61.9 cm (24.4")
SEATING SIZES (TRU-BALANCE® 3 and 4)	
TRU-Balance®	Width: 40.5cm-55.8cm (16-22") Depth: 30.5cm-55.8cm (16-22")
Q Captain TRU-Balance® 3 Static Seat	Width: 30.5cm-55.8cm (12-22") Depth: 25.4cm-55.8 cm (10-22")
TRU-Balance® 3 Power Tilt	Width: 30.5cm-55.8cm (12-22") Depth: 25.4cm-55.8cm (10-22")
TRU-Balance® 3 Power Lift & Tilt	Width: 30.5cm-55.8cm (12-22") Depth: 25.4cm-55.8cm (12-22")
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance®4	44.5cm - 48.9cm (17.5 -19.25")
Q Captain / Static	41.9cm - 47cm (16.5 -18.5")
TRU-Balance® 3 Power Tilt ^{5/7}	44.5cm - 49.5cm (17.5 -19.5")
TRU-Balance® 3 Tilt & iLevel® ^{5/7}	44.5 cm - 49.5cm (17.5"- 19.5")
BATTERY SIZE	MK 50Ah(2)
BATTERY WEIGHT ⁶	17.2 kg each
AVAILABLE ELECTRONICS	120A Q-Logic 3
RANGE PER CHARGE ³	Up to 32 km
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-pole or 4-poles motors
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT	193 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 136 kg. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 1.5kg. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded.

Standard Features

- **Upgraded** Independent SRS (Smooth Ride Suspension)
- Fender lights & USB Charger
- 7.24 km/h at i level®
- Matching wheel inserts
- Made in Australia of foreign and domestic parts
- Available with optional Accu-Trac Advanced Tracking Technology

Options

- 12" iLevel® Power seat elevation allows driving at speeds up to 7.2 km/h while elevated

Colors

