

# Q6 EDGE 3 *Stretto*

Q6 EDGE 3 STRETTO™  
shown with iLevel® power seat elevation,  
Q-Logic 3 environmental controls, and  
Stealth® TRU-Comfort 2 cushion.



## BUILT COMPACT

The Q6 Edge® 3 Stretto™ features an ultra-narrow base for exceptional maneuverability in tight indoor spaces. It is equipped with independent SRS (Smooth Ride Suspension), making it an excellent choice for individuals who need a compact yet capable Mid-Wheel Drive. It's available with 31.8 cm or 35.6 cm drive wheels, allowing for a customised fit, and offers the option of 30.5 cm (12") of iLevel® power adjustable seat height for improved access, interaction, and independence.



| [quantumrehab.com.au](http://quantumrehab.com.au)

## Specifications

<b>DRIVE WHEELS (Mid-Wheel Drive)</b>	35.5 cm (14") or 31.7 cm (12.5") Pneumatic or Flat Free Tyres
<b>CASTER WHEELS</b>	<b>Front:</b> 12.7 cm or 15.2cm (5") Solid <b>Rear:</b> 15.2 cm (6") Solid
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 10 km/h Up to 5.6 km/h w/ iLevel activated
<b>GROUND CLEARANCE</b>	7.3 cm
<b>TURNING RADIUS</b>	53.3 cm
<b>OVERALL LENGTH <sup>2</sup></b>	82.5 cm without front riggings
<b>BASE WIDTH</b>	54.6 cm with 35.5 cm drive wheels 51.9 cm with 31.5 cm drive wheels
<b>SEATING SIZES ( TRU-BALANCE® 3 and 4 )</b>	
<b>TRU-Balance 4 ®</b>	<b>Width:</b> 40.5 cm - 50.8 cm (16-20") <b>Depth:</b> 40.5 cm - 50.8 cm (16-20")
<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 30.5 cm - 50.8cm (12-20") <b>Depth:</b> 30.5cm - 50.8cm (12-20")
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 30.5 cm - 50.8cm (12-20") <b>Depth:</b> 30.5 cm - 50.8cm (12-20")
<b>TRU-Balance® 3 Power Lift &amp; Tilt</b>	<b>Width:</b> 30.5 cm - 50.8cm (12-20") <b>Depth:</b> 30.5 cm - 50.8cm (12-20")
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>TRU-Balance®4</b>	18-19" (45.7cm - 48.2cm)
<b>TRU-Balance® 3 Static</b>	16-18" (40.5cm - 45.7cm)
<b>TRU-Balance® 3 Power Tilt <sup>5/7</sup></b>	17.5-19.5" (44.45cm - 49.5cm)
<b>TRU-Balance® 3 Tilt &amp; iLevel® <sup>5/7</sup></b>	17.5-19.5" (44.5cm - 49.5cm)
<b>BATTERY SIZE</b>	
<b>BATTERY WEIGHT <sup>6</sup></b>	MK 48Ah(2) 17 kg each
<b>AVAILABLE ELECTRONICS</b>	120A Q-Logic 3
<b>RANGE PER CHARGE <sup>3</sup></b>	Up to 25.6 km
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-pole motors
<b>WEIGHT CAPACITY <sup>4</sup></b>	136 kg
<b>BASE WEIGHT</b>	185 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 136 kg. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.  
Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 1.5kg. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded.
- 10) Limits front rigging options


## Standard Features

- 2-pole motors
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Q-Logic 3 Drive controls
- Speeds up to 10 kph
- LED front and rear fender lights
- Quiet ride technology
- USB charger included
- ISO 7176-19 transit crash tested approved with all configurations

## Options

- TRU-Balance® 4 Power positioning:  
10 degree or 30 degree anterior tilt Options with memory seating
- TRU-Balance® 3 Power Positioning:  
iLevel® option allows 30.5 cm 12" of seat elevation while at walking speed (5.6 km/h)
- Power Tilt, Recline, and Elevating Legrests that can power to Floor for easy transfers.
- Accu -Trac advanced tracking technology
- Optional 6" front castors<sup>10</sup>
- LED Lighting Package
- NEW backup camera

## Colors



**QUANTUM®**  
#1 FOR REHAB POWER

 **PROUDLY BUILT IN AUSTRALIA**

(03) 8770 9600 | [quantumrehab.com.au](http://quantumrehab.com.au)